

MODEL NUMBER	INTERNAL DIMS.		ZONE LENGTH			EXTERNAL DIMS.			POWER		ZONES WATTS		
	ID	L	1	2	3	H	W	L	TOTAL		1	2	3
									WATTS	VOLTS			
18" HEATED LENGTH													
SST- 0.75 -0- 18 -3C- D-2155- R	0.750 x 18		3	12	3	14.5 x 10.5 x 24			855	115	142.5	570	142.5
SST- 1.00 -0- 18 -3C- D-2155- S	1.000 x 18		3	12	3	14.5 x 10.5 x 24			1140	115	190	760	190.0
SST- 1.25 -0- 18 -3C- D-2155- T	1.250 x 18		3	12	3	14.5 x 10.5 x 24			1800	230	300	1200	300.0
SST- 1.63 -0- 18 -3C- D-2155- U	1.625 x 18		3	12	3	14.5 x 10.5 x 24			2160	230	360	1440	360.0
SST- 2.00 -0- 18 -3C- D-2155- V	2.000 x 18		3	12	3	17.0 x 13.0 x 24			2475	230	412.5	1650	412.5
SST- 2.38 -0- 18 -3C- D-2155- W	2.375 x 18		3	12	3	17.0 x 13.0 x 24			2850	230	475	1900	475.0
SST- 3.00 -0- 18 -3C- D-2155- Y	3.000 x 18		3	12	3	17.0 x 13.0 x 24			3390	230	565	2260	565.0
SST- 3.75 -0- 18 -3C- D-2155- CA	3.750 x 18		3	12	3	20.5 x 16.5 x 24			4950	230	825	3300	825.0
SST- 4.13 -0- 18 -3C- D-2155- CB	4.125 x 18		3	12	3	20.5 x 16.5 x 24			5100	230	850	3400	850.0
SST- 5.00 -0- 18 -3C- D-2155- CC	5.000 x 18		3	12	3	20.5 x 16.5 x 24			6030	230	1005	4020	1005
SST- 5.50 -0- 18 -3C- D-2155- CD	5.500 x 18		3	12	3	20.5 x 16.5 x 24			6690	230	1115	4460	1115
SST- 6.00 -0- 18 -3C- D-2155- CE	6.000 x 18		3	12	3	20.5 x 16.5 x 24			6900	230	1150	4600	1150

24" HEATED LENGTH

SST- 0.75 -0- 24 -3C- D-2155- AG	0.750 x 24		6	12	6	14.5 x 13.0 x 30			1140	115	285	570	285
SST- 1.00 -0- 24 -3C- D-2155- AH	1.000 x 24		6	12	6	14.5 x 13.0 x 30			1520	230	380	760	380
SST- 1.25 -0- 24 -3C- D-2155- AJ	1.250 x 24		6	12	6	14.5 x 13.0 x 30			2400	230	600	1200	600
SST- 1.63 -0- 24 -3C- D-2155- AK	1.625 x 24		6	12	6	14.5 x 13.0 x 30			2880	230	720	1440	720
SST- 2.00 -0- 24 -3C- D-2155- AL	2.000 x 24		6	12	6	17.0 x 13.0 x 30			3300	230	825	1650	825
SST- 2.38 -0- 24 -3C- D-2155- AM	2.375 x 24		6	12	6	17.0 x 13.0 x 30			3800	230	950	1900	950
SST- 3.00 -0- 24 -3C- D-2155- AN	3.000 x 24		6	12	6	17.0 x 13.0 x 30			4520	230	1130	2260	1130
SST- 3.75 -0- 24 -3C- D-2155- CL	3.750 x 24		6	12	6	22.5 x 16.5 x 30			6600	230	1650	3300	1650
SST- 4.13 -0- 24 -3C- D-2155- CM	4.125 x 24		6	12	6	22.5 x 16.5 x 30			6800	230	1700	3400	1700
SST- 5.00 -0- 24 -3C- D-2155- CN	5.000 x 24		6	12	6	22.5 x 16.5 x 30			8040	230	2010	4020	2010
SST- 5.50 -0- 24 -3C- D-2155- CP	5.500 x 24		6	12	6	22.5 x 16.5 x 30			8960	230	2250	4460	2250
SST- 6.00 -0- 24 -3C- D-2155- CR	6.000 x 24		6	12	6	22.5 x 16.5 x 30			9200	230	2300	4600	2300

30" HEATED LENGTH

SST- 0.75 -0- 30 -3C- D-2155- AY	0.750 x 30		6	18	6	14.5 x 13.0 x 36			1430	230	290	850	290
SST- 1.00 -0- 30 -3C- D-2155- AZ	1.000 x 30		6	18	6	14.5 x 13.0 x 36			1900	230	380	1140	380
SST- 1.25 -0- 30 -3C- D-2155- BA	1.250 x 30		6	18	6	14.5 x 13.0 x 36			2418	230	484	1450	484
SST- 1.63 -0- 30 -3C- D-2155- BB	1.625 x 30		6	18	6	14.5 x 13.0 x 36			3366	230	673	2020	673
SST- 2.00 -0- 30 -3C- D-2155- BC	2.000 x 30		6	18	6	17.0 x 13.0 x 36			3984	230	797	2390	797
SST- 2.38 -0- 30 -3C- D-2155- BD	2.375 x 30		6	18	6	17.0 x 13.0 x 36			4316	230	863	2590	863
SST- 3.00 -0- 30 -3C- D-2155- BE	3.000 x 30		6	18	6	17.0 x 13.0 x 36			5734	230	1147	3440	1147
SST- 3.75 -0- 30 -3C- D-2155- CY	3.750 x 30		6	18	6	22.5 x 16.5 x 36			7432	230	1486	4460	1486
SST- 4.13 -0- 30 -3C- D-2155- CZ	4.125 x 30		6	18	6	22.5 x 16.5 x 36			7500	230	1500	4500	1500
SST- 5.00 -0- 30 -3C- D-2155- DA	5.000 x 30		6	18	6	22.5 x 16.5 x 36			7368	230	1674	4020	1674
SST- 5.50 -0- 30 -3C- D-2155- DB	5.500 x 30		6	18	6	22.5 x 16.5 x 36			8960	230	2250	4460	2250
SST- 6.00 -0- 30 -3C- D-2155- DC	6.000 x 30		6	18	6	22.5 x 16.5 x 36			9534	230	2467	4600	2467

36" HEATED LENGTH

SST- 0.75 -0- 36 -3C- D-2155- BN	0.750 x 36		6	12-12	6	14.5 x 13.0 x 42			1710	230	285	1140	285
SST- 1.00 -0- 36 -3C- D-2155- BP	1.000 x 36		6	12-12	6	14.5 x 13.0 x 42			2280	230	380	1520	380
SST- 1.25 -0- 36 -3C- D-2155- BR	1.250 x 36		6	12-12	6	14.5 x 13.0 x 42			3600	230	600	2400	600
SST- 1.63 -0- 36 -3C- D-2155- BS	1.625 x 36		6	12-12	6	14.5 x 13.0 x 42			4320	230	720	2880	720
SST- 2.00 -0- 36 -3C- D-2155- BT	2.000 x 36		6	12-12	6	17.0 x 13.0 x 42			4950	230	825	3300	825
SST- 2.38 -0- 36 -3C- D-2155- BU	2.375 x 36		6	12-12	6	17.0 x 13.0 x 42			5700	230	950	3800	950
SST- 3.00 -0- 36 -3C- D-2155- BV	3.000 x 36		6	12-12	6	17.0 x 13.0 x 42			6780	230	1130	4520	1130
SST- 3.75 -0- 36 -3C- D-2155- DJ	3.750 x 36		6	12-12	6	22.5 x 16.5 x 42			9900	230	1650	6600	1650
SST- 4.13 -0- 36 -3C- D-2155- DK	4.125 x 36		6	12-12	6	22.5 x 16.5 x 42			10200	230	1700	6800	1700
SST- 5.00 -0- 36 -3C- D-2155- DL	5.000 x 36		6	12-12	6	22.5 x 16.5 x 42			12060	230	2010	8040	2010
SST- 5.50 -0- 36 -3C- D-2155- DM	5.500 x 36		6	12-12	6	22.5 x 16.5 x 42			13420	230	2250	8920	2250
SST- 6.00 -0- 36 -3C- D-2155- DN	6.000 x 36		6	12-12	6	22.5 x 16.5 x 42			13800	230	2300	9200	2300